#### 1 Introduction

This document describes the Printheader class, which wraps the program printheader to provide information about the image headers. This will be used in the Dataset object.

## 2 Input

All this takes as input is a frame set in setImage().

### 3 Output

The following results are produced:

```
{'exposure_time': 5.0,
'distance': 170.0,
'phi_start': 290.0,
'phi_end': 291.0,
'epoch': 1096203695.0,
'date': 'Sun Sep 26 14:01:35 2004',
'wavelength': 0.979660,
'detector': 'adsc'}
```

# 4 Updates

## 4.1 15/JUN/06

Added a HeaderCache singleton to record the image headers as they are read, so that when they are asked for again the results come back much more rapidly.