

1 Introduction

The application `xia2scan` is designed to scan through a list of provided images and print information about the “quality” of the diffraction on the images. This will end up as an executable `xia2-scan`(=> `.sh/.bat`) and a python program `xia2scan.py`.

1.1 Appication

This is designed to be useful for

- Selecting the best images for autoindexing from a sweep.
- Optimizing humidity when a humidifier is available.

by printing information about the strength of diffraction in each image.

2 Uses...

This uses the wrappers for `labelit.screen` and `labelit.stats_distl`.

3 Dependencies

This application will depend on having access to a user provided beam centre (by implication then an input argument / command line handler) and also diffraction image name parsing to allow searching of a directory for matching images.

4 Use Cases

4.1 UC 1: Humidifier

Requires:

- The correct beam position (for indexing from single images.)
- The images.

Provides:

- The statistics for each image in the list.

4.2 UC 2: Indexing Image Selection

This will scan through all of the images in a sweep and select, based on the spot separation and so on, the best pair for autoindexing.

5 Implementation

5.1 UC 1: Humidifier

This will use the following classes:

- Schema.Sweep.SweepFactory - to generate a list of sweeps from the images provided.
- Handlers.CommandLine - to hold the beam centre information from a previous indexing run.

To achieve reasonable speed, the labelit run will not optimise the beam centre (this is unreliable from a single image anyway) and also will run the labelit before distl, and use labelit.stats.distl to get the results afterwards.

This should do all of the running, then collate the results and print a summary. The following arc should be followed:

- [determine accurate beam centre]
- Digest matching frames to list of sweeps.
- For each sweep, run labelit, get an idea of the unit cell volume, run labelit.stats.distl to get the intensity summary, save.
- Print summary by image number.

5.1.1 Test Data

Rajan provided me with a couple of test data sets - use these. Oh bugger - the images are stills... Add a special case in the sweep factory to return these as separate frames? Or just plough ahead??